


# Discover Faster Publish Quicker with Moxi GO™

Run 15 assays in under 15 minutes with simultaneous cell counts


- Configurable 4-Plex Bead ELISA's
- CRISPR/Transfection checks
- In cell Western's
- Rapid phenotyping, immuno-profiling
- Gold standard cell count and viability
- No training, maintenance, or cleaning



**Specs:** 1 laser (532/488) • 1 PMT (FITC/PE) • Count and Cell Volume (Coulter Principle) • FCS/Bitmap/CSV file output  
USB-On-The-Go • Dual Use Flow Cell = No cleaning, maintenance or cross contamination • 10k-IM cells/ml




1. Run Control



2. Run Treated Samples

Bead	Plex-002		Ratio
	Reference MFI (mv)	Unknown MFI (mv)	
1	13	1423	113.40
2	31	976	31.84
3	36	1898	52.37
4	94	89	0.94

3. Summary Plex Table



4. Relative Quant Readout

Product	Qty.	ORFLO Cat. No.
Moxi GO Kit	1	MXG001
<b>Accessories</b>		
Moxi Flow Cassettes, Type MF-S	Pack	MXC020
Moxi Flow Cassettes, Type MF-M	Pack	MXC021
Moxi Cyte Viability Reagent	Bottle (15ml, 100 tests)	MXA055
Moxi Flow System Check Beads, 5mL	Bottle (5ml)	MXA016
Moxi Cyte Flow Reagent	Bottle (15ml, 750 tests)	MXA080

Purchase before  
2017 and receive  
\$3,000 off the  
(\$12,995 - \$14,995)  
MSRP

© Copyright 2016 Orflo Technologies. All rights reserved.

www.orflo.com

855.879.6694

Next Gen  
Flow  
Cytometry™

